

FETAL ECG DETECTION METHOD FOR ABDOMINAL RECORDINGS

**final project presentation for
Biomedical Signal Processing**

Zofia Mizgalewicz

08.02.2024



SELECTED ARTICLE

IOP PUBLISHING

PHYSIOLOGICAL MEASUREMENT

Physiol. Meas. **28** (2007) 373–388

[doi:10.1088/0967-3334/28/4/004](https://doi.org/10.1088/0967-3334/28/4/004)

A robust fetal ECG detection method for abdominal recordings

**Suzanna M M Martens¹, Chiara Rabotti¹, Massimo Mischi¹ and
Rob J Sluijter²**

¹ Department of Electrical Engineering, University of Technology Eindhoven, Eindhoven,
The Netherlands

² Philips Research Laboratories, Eindhoven, The Netherlands

E-mail: smm.martens@gmail.com

Received 29 November 2006, accepted for publication 8 February 2007

Published 7 March 2007

Online at stacks.iop.org/PM/28/373

DATASET

- 25 recordings of abdominal ECG
- 4 channel ECG
- each 1 minute long
- sampling frequency of 1000Hz



PhysioNet

© Challenge  Open Access

Noninvasive Fetal ECG: The PhysioNet/Computing in Cardiology Challenge 2013

Ikaro Silva , Joachim Behar , Reza Sameni , Tingting Zhu , Gari D. Clifford , George Moody 

PIPELINE

01

DATA LOADING AND PREPROCESSING

02

QRS DETECTOR

03

MECG CANCELLATION METHOD

04

FECG FEATURE EXTRACTION

05

EVALUATION

06

RESULTS





Figure 1. Nomenclature of the electrocardiogram (ECG).

DATA PREPROCESSING

01

LOADING AND
SIGNAL
INTERPOLATION

02

BASELINE
WANDERING
REMOVAL

high pass FIR filter
cutoff 3Hz
1001 taps
sampling frequency 1000Hz

03

POWER-LINE
INTERFERENC
E CANCELLER


notch filter
cutoff 50Hz
fs 1000Hz

04

UPSAMPLING
OF THE
SIGNAL

from 1000 Hz
to 2000 Hz


FILTERING




Baseline wander removal. Left plot shows part of one abdominal ECG channel before and after baseline wander removal with the filter. The right plot shows the amplitude (left) and phase (right) of the high-pass filter transfer function.

VISUAL INSPECTION

Original signal for patient a13



Filtered signal for patient a13



comparison of the original signal and signal after preprocessing

QRS DETECTOR

optimal combination of the channels is derived by performing a principal component analysis

ENHANCE THE QRS SIGNAL

minimal distance between R peaks:



- 0.6s for MECCG (HR 100bpm)
- 0.34s for FECCG (FHR 175bpm)

take a bunch of segments of the enhanced signal that contains a QRS complex and take the mean that will serve as a QRS template

CREATE A QRS TEMPLATE


use cross-correlation between the enhanced signal and the QRS template to identify potential QRS complexes

DETECT QRS PEAKS





QRS detection for the MECG (left plot) and FECG (right plot). Plots show signals from 4 leads and the signal generated by the multi-channel QRS enhancement method (bottom signal). The detected QRS complexes are indicated by purple x.

MECG CANCELLER




FECG DETECTOR


averages FECG complexes synchronized on the fetal QRS complex
returns FECG complex template and FHR trace



EVALUATION



RESULTS



```
Accuracy for patient 1: 56.72%
Accuracy for patient 2: 50.77%
Accuracy for patient 3: 82.81%
Accuracy for patient 4: 95.45%
Accuracy for patient 5: 100.00%
Accuracy for patient 6: 49.21%
Accuracy for patient 7: 56.35%
Accuracy for patient 8: 97.66%
Accuracy for patient 9: 54.47%
Accuracy for patient 10: 56.91%
Accuracy for patient 11: 58.73%
Accuracy for patient 12: 84.78%
Accuracy for patient 13: 99.21%
Accuracy for patient 14: 80.31%
Accuracy for patient 15: 75.37%
Accuracy for patient 16: 49.12%
Accuracy for patient 17: 96.97%
Accuracy for patient 18: 50.00%
Accuracy for patient 19: 94.49%
Accuracy for patient 20: 95.45%
Accuracy for patient 21: 66.91%
Accuracy for patient 22: 60.28%
Accuracy for patient 23: 98.43%
Accuracy for patient 24: 99.19%
Accuracy for patient 25: 98.40%
Mean accuracy: 76.32%
```